


# Olecranon Hook Plate™


Elbow & Forearm System






### Exposure and Reduction

- Expose the olecranon with a standard posterior approach.
- Reduce the fracture and hold with K-wire(s).
- Make two longitudinal incisions in line with the fibers of the triceps insertion down to the proximal tip of the olecranon.



### Guide Preparation

- Apply the Hook Drill Guide to the olecranon ensuring the guide is completely seated and flush with the bone. Drill two holes through the outer cortex with a 2.3mm (red) drill.




### Hook Plate Insertion

- Insert the plate's hooks into the holes drilled in the olecranon.
- Seat the Impactor between the two hooks and impact the plate into the distal fragment.
- Confirm that the reduction is adequate.
- Place a non-locking 3.2mm cortical screw in the end of the slotted hole that is farthest from the fracture.



### Expander/Compression Tool

- Insert screwdriver tip of the Expander/Compression Tool into screw head and engage jaw into the adjacent screw hole away from the fracture or osteotomy.
- Loosen screw **1/4** turn.
- Gently squeeze the Expander/Compression Tool to obtain the desired compression. Tighten screw.



### Final Fixation

- Complete fixation with 3.2mm cortical screws, as needed.

All implants made from surgical grade stainless steel

## Olecranon Hook Plate™

*Neutral*  
OH00K-4


*Left*  
OH00KL-6  
OH00KL-8

*Right*  
OH00KR-6  
OH00KR-8



## Expander/ Compression Tool

XPANDR



## Olecranon Hook Plate™ Drill Guide

GDOH00K-2.3



## Drill Guide

GUIDELCBS-2.3



## Cortical Screw

HEX3.2-xx  
08mm to 40mm



## Cortical Locking Screw

LCBS3.2-xx  
08mm to 30mm



## Impactor

IMPCT-SLED



TriMed, Inc. / 27533 Avenue Hopkins / Valencia, CA 91355 USA / 800-633-7221 / [www.trimedortho.com](http://www.trimedortho.com)

The technique presented is one suggested surgical technique. The decision to use a specific implant and the surgical technique must be based on sound medical judgment by the surgeon that takes into consideration factors such as the circumstances and configuration of the injury.

This document is controlled by TriMed, Inc. When downloaded, printed, and/or copied, this document becomes uncontrolled, and users should always check [trimedortho.com](http://trimedortho.com) for the latest version.

For indications, contraindications, warnings and precautions related to TriMed Elbow & Forearm System reference IFU on [trimedortho.com/ifu](http://trimedortho.com/ifu).

See [trimedortho.com/patents](http://trimedortho.com/patents) for all patent information.