



# Total Wrist Fusion Plate

Total Wrist Fusion System




SURGICAL TECHNIQUE




## Preparation

- Expose and decorticate the articular surfaces for the size and shape of plate to be used.
- If necessary, create a trough to allow plate to sit flush with the bone.




## Distal Fixation (2.7mm Screws)

- Complete all fixation distal to the radiocarpal joint.
- Use 2.0mm (yellow) drill for 2.7mm screws.
  - Use GUIDELCBS-2.0 for locking screws.
  - Use GUIDE-2.0/2.7 for non-locking screws



## Technique for Adding Compression

- Place a 3.2mm cortical screw in the end of the slotted hole farthest from the joint. Do not completely seat against the plate.
- Engage the driver tip of the Expander/Compression Tool into the screw head. Insert the hook of the Expander/Compression Tool in the adjacent empty screw hole away from the joint.
- Gently squeeze the Expander/Compression Tool to obtain desired compression. Tighten the screw.



## Proximal Fixation (3.2mm Screws)

- Complete proximal fixation in the normal manner.
- Use 2.3mm (red) drill for 3.2mm screws.
  - Use GUIDELCBS-2.3 for locking screws.
  - Use GUIDE-2.3/3.2 for non-locking screws

TriMed, Inc. / 27533 Avenue Hopkins / Santa Clarita, CA 91355 USA / 800-633-7221 / [www.trimedortho.com](http://www.trimedortho.com)

The technique presented is one suggested surgical technique. The decision to use a specific implant and the surgical technique must be based on sound medical judgment by the surgeon that takes into consideration factors such as the circumstances and configuration of the injury.

This document is controlled by TriMed, Inc. When downloaded, printed, and/or copied, this document becomes uncontrolled, and users should always check [trimedortho.com](http://trimedortho.com) for the latest version.

For indications, contraindications, warnings and precautions related to TriMed Total Wrist Fusion System reference IFU on [trimedortho.com/ifu](http://trimedortho.com/ifu).

See [trimedortho.com/patents](http://trimedortho.com/patents) for all patent information.